VAPOTHERM®

# Nitric Oxide Disposable Patient Circuit

## Instructions For Use

For further information contact:

Vapotherm Inc.  
100 Domain Drive  
Exeter, NH 03833  
USA  
Phone: 603-658-0011  
Fax: 603-658-0181  
www.vapotherm.com

May be patented  
[www.vapotherm.com/patents](http://www.vapotherm.com/patents)

Technical Support Line  
Domestic: 855-557-8276  
International: 603-658-5121  
TS@Vtherm.com

RMS – UK Limited  
28 Trinity Road  
Nailsea, North Somerset BS48 4NU  
United Kingdom  
Phone: +44-1275-85-88-91  
Fax: +44-1275-85-88-91

Vapotherm Inc. has declared that this product conforms with the European Council Directive 93/42/EEC Medical Device Directive when it is used in accordance with the instructions provided in the Instructions For Use.

## Precision Flow® Disposable Patient Circuit for use with the Nitric Oxide Delivery System

### Nitric Oxide Delivery System

#### Operations with the Nitric Oxide Delivery System

1. Assemble and then install the Disposable Patient Circuit (DPC) in the Precision Flow unit in accordance with Vapotherm instructions.
2. Attach the 22mm and 15mm tubing connectors of the DPC to the Injector Module and let hang below the Precision Flow unit as shown in **Figure 1**.
3. Verify the tubing connectors are pressed firmly into the Injector Module to prevent leaks.
4. Connect the sample tee per instructions. Please refer to your Nitric Oxide Delivery System Operator's Manual.
5. Route the tubing so it exits the Precision Flow at the bottom of the DPC next to the delivery tube so that it does not interfere with the door closing (**Figure 1**).
6. Close the Precision Flow door, making sure the tubing does not interfere, kink, or bind in the door (**Figure 2**). This may stop the flow through the circuit and result in an occluded tube alarm.
7. Connect cannula to the end of the gas sample tee. Verify flow is passing through the cannula before attaching it to the patient.
8. Allow the Precision Flow to heat to set temperature before placing cannula on the patient or administering Inhaled Nitric Oxide (INO) therapy.
9. Follow the Nitric Oxide Delivery System Information for Use for proper operation of the Nitric Oxide Delivery System equipment and nitric oxide delivery.


The Nitric Oxide Delivery System adds NO/N<sub>2</sub> gas flow to the breathing circuit flow in proportion to the NO setting (up to 10% at 80ppm), and subtracts gas from the breathing circuit via the gas sampling system at a nominal flow rate of 0.23 L/min. This slightly reduces the delivered gas to the patient - you may wish to adjust the gas source flow rate to compensate.

#### High Flow Therapy after discontinuing Nitric Oxide


Remove the tubing connectors from the Injector Module and attach the two connectors together (**see Figure 3**). The connectors are compatible and will complete the air path to permit operation without Nitric Oxide Delivery System equipment. Flow will exit the injector port opening of the Injector Module if left in series with the circuit. Remove the cannula from the gas sample tee, remove the tee from the system and connect the cannula to the end of the DPC delivery tube.

**Warning:** Verify the connectors are firmly pressed together and the cannula is firmly connected to the DPC delivery tube. Verify that there are no leaks in the circuit. If the connections are not properly attached, flow through the DPC may be inaccurate or non-existent.


The Precision Flow® has been validated for use with the INOmax and NOxBox Systems. Please refer to your Nitric Oxide Delivery System Operator's Manual.

**Figure 1**  
Door Open

**Figure 2**  
Door Closed

**Figure 3**  
Removing Injector Module

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