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43000719-FR-FR Rev D

Aerosol Adapter - Instructions


Intended Use:
This Vapotherm Aerosol Adapter AAA-2 is intended to facilitate the connection between Vapotherm High Velocity Therapy HVT® 2.0 and the Aerogen Solo Nebulizer.

WARNINGS:

- The AAA-2 is a single-patient use product that can be used for up to 30 days from the date of opening the packaging, replace if the adapter appears soiled.
- Do Not use the device on patient for longer than 30 days.
- Clinicians should follow hospital best practice for infection control and AAA-2 removal storage & labeling for intermittent use.
- It is important to maintain proper upright orientation of the patient adapter and nebulizer unit during the drug administration process to prevent rainout from reaching the patient.
- Lower flows with high velocity therapy may require more frequent condensate drainage.
- Failure to drain condensate may result in condensate delivery to the patient.
- The adapter must be removed between aerosol treatments, as condensate will accumulate if left in the circuit for extended periods.
- Ensure that the patient interface and delivery tube are reconnected whenever the adapter is removed to continue with high velocity therapy.

Instructions for Use:

1. Use an aseptic technique to place the (1) Aerosol Adapter in line between the patient delivery tube and the (2) ProSoft cannula or the Vapotherm TA-22 (tracheostomy tube adapter).
2. Set the HVT 2.0 device to a recommended flow rate in the range of 10-20 Liters per Minute.
3. Disconnect the (2) patient interface from the (3) delivery tube.
4. Gently connect the adapter to the delivery tube to facilitate easy removal.
5. Then, connect the (2) patient interface to the (1) adapter.
6. Once the (1) adapter is in place, connect the (4) nebulizer components and start the aerosol treatment. **Note:** It is important to **maintain a proper upright orientation** of the (1) patient adapter and (4) nebulizer unit during the drug administration process.
7. After the nebulization treatment is complete, remove the (1) adapter from (2) patient interface and the (3) delivery tube first by **holding the (5) rubber bushing on the (1) adapter** and then twisting/pulling the (3) delivery tube out. Then remove the (2) patient interface from the (1) adapter using the same twisting/pulling motion.

- (1) Aerosol Adapter
- (2) Patient Interface (ProSoft 2.0 Cannula or TA-22)
- (3) Delivery Tube
- (4) Aerogen Solo Nebulizer
- (5) Rubber Bushing

Note: If the (5) rubber bushing comes off when disconnecting the adapter from the delivery tube, manually push the bushing back into the adapter as seen in the image above.

8. Reconnect the (2) patient interface to the (3) delivery tube.
9. **Empty or drain** any condensate collected in the adapter per your institution's policy.

Refer to the Aerogen SOLO Instructions for Use for aerosol nebulization instructions.

Refer to the HVT® 2.0 Instructions for Use for the operation of the HVT® 2.0.

43000719-EN-US Rev D

VAPOTHERM® Aerosol Adapter – Instructions


Intended Use:
This Vapotherm Aerosol Adapter AAA-2 is intended to facilitate the connection between Vapotherm High Velocity Therapy HVT® 2.0 and the Aerogen Solo Nebulizer.

WARNINGS:

- The AAA-2 is a single-patient use product that can be used for up to 30 days from the date of opening the packaging, replace if the adapter appears soiled.
- Do Not use the device on patient for longer than 30 days.
- Clinicians should follow hospital best practice for infection control and AAA-2 removal storage & labelling for intermittent use.
- It is important to maintain proper upright orientation of the patient adapter and nebulizer unit during the drug administration process to prevent rainout from reaching the patient.
- Lower flows with high velocity therapy may require more frequent condensate drainage.
- Failure to drain condensate may result in condensate delivery to the patient.
- The adapter must be removed between aerosol treatments, as condensate will accumulate if left in the circuit for extended periods.
- Ensure that the patient interface and delivery tube are reconnected whenever the adapter is removed to continue with high velocity therapy.

Instructions for Use:

- Use an aseptic technique to place the (1) Aerosol Adapter in line between the patient delivery tube and the (2) ProSoft cannula or the Vapotherm TA-22 (tracheostomy tube adapter).
- Set the HVT 2.0 device to a recommended flow rate in the range of 10-20 Litres per Minute
- Disconnect the (2) patient interface from the (3) delivery tube.
- Gently connect the adapter to the delivery tube to facilitate easy removal.
- Then, connect the (2) patient interface to the (1) adapter.
- Once the (1) adapter is in place, connect the (4) nebuliser components and start the aerosol treatment. **Note:** It is important to maintain a proper upright orientation of the (1) patient adapter and (4) nebuliser unit during the drug administration process.
- After the nebulisation treatment is complete, remove the (1) adapter from (2) patient interface and the (3) delivery tube first by holding the (5) rubber bushing on the (1) adapter and then twisting/pulling the (3) delivery tube out. Then remove the (2) patient interface from the (1) adapter using the same twisting/pulling motion.

- (1) Aerosol Adapter
- (2) Patient Interface (ProSoft 2.0 Cannula or TA-22)
- (3) Delivery Tube
- (4) Aerogen Solo Nebulizer
- (5) Rubber Bushing

Note: If the (5) rubber bushing comes off when disconnecting the adapter from the delivery tube, manually push the bushing back into the adapter as seen in the image above.

- Reconnect the (2) patient interface to the (3) delivery tube.
- Empty or drain any condensate collected in the adapter per your institution's policy.

Refer to the Aerogen SOLO Instructions for Use for aerosol nebulisation instructions.

Refer to the HVT® 2.0 Instructions for Use for the operation of the HVT® 2.0.

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