


<b>Technical Bulletin</b>		
 <b>VAPOTHERM</b>	<b>TITLE: HVT 2.0 Instructions for Use Updates</b>	
DOC #: TECH-00045_EN-GB	Revision: A	PAGE: 1 of 1

**Distribution:** Vapotherm website, Customer Service, Technical Support and Sales Team (Global)

**Subject:** This bulletin serves as an update to your HVT 2.0 Instructions for Use:

- 43000655 Rev. B (EN)
- 43000656 Rev. B (EN)
- 43000657 Rev. A (DE)
- 43000671 Rev. A (GB)
- 43000656 Rev. A (FR, IT, PT-BR, NL-BE, SV, PT, TR, ES-MX, ES, ID, JA-JP)

See the referenced Instructions for Use (IFU) listed above for complete instructions and all warnings and precautions. IFU's can be downloaded at <https://vapotherm.com/hvt2/>

**Date:** Effective Immediately

Please remove or replace the following information within the Sections of the IFU as indicated:

### Section 4: Using the HVT 2.0 System

The Pediatric Flow Calculator feature is no longer available on the HVT 2.0 device.

### Appendix C – Technical Specifications

#### System Requirements

AC power:	100 to 240VAC, 50 to 60Hz
Oxygen supply:	Oxygen or oxygen concentration at inlet pressures between 4 and 87 psi.

#### Performance

Oxygen percentage:	Accuracy $\pm$ 2% Wall or Tank Source
--------------------	---------------------------------------

If you have any questions about this change, please contact Vapotherm Technical Support at [ts@vtherm.com](mailto:ts@vtherm.com) or get in touch with your local Vapotherm representative.

AJW Technology Consulting GmbH  
Breite Straße 3  
40213 Düsseldorf (Germany)  
Telephone: +49 211 54059 6030