


<b>Technical Bulletin</b>		
 <b>VAPOTHERM</b>	<b>TITLE: HVT 2.0 Instructions for Use Updates</b>	
DOC #: TECH-00048	Revision: A	PAGE: 1 of 1

**Distribution:** Vapotherm website, Customer Service, Technical Support, and Sales Team (Global)

**Subject:** This Bulletin serves as an update to your HVT 2.0 Instructions for Use:

- 43000655 Rev. B (US)
- 43000656 Rev. B (International)

See the referenced Instructions for Use (IFU) listed above for complete instructions and all warnings and precautions. IFUs can be downloaded at <https://vapotherm.com/hvt2/>

**Date:** January 25, 2023

Please replace the following information within the Section(s) of the IFU as indicated:

## Appendix A: Nurse Call System Installations

### Nurse Call Hardware Interface Description

The Nurse Call Communication Cable is available in four variants:

Cable Type
NO contact
NO + 10k contact
NC contact
NC+10k contact

Please order the cable type compatible with your nurse call system.

If you have any questions about this change, please contact Vapotherm Technical Support at [ts@vtherm.com](mailto:ts@vtherm.com) or your local Vapotherm representative.